




Cleaning the Spiriva HandiHaler



- Clean the HandiHaler once a month or as needed
 - Open the lid and mouthpiece
 - Open the base by lifting the green piercing button
 - Rinse the HandiHaler with warm water to remove any powder
- Do not use soap or place in the dishwasher, just rinse with water
- Dry the HandiHaler completely, tip out the excess water and leave the lid, mouthpiece, and base open
- The best time to clean the HandiHaler is right after use
- The HandiHaler takes 24 hours to air dry

May 2007